Ophthalmology Times Research Scholar Honoree Program

Application of Active Learning of the Contrast Sensitivity Function in Patients with Macular Disease and Good Visual Acuity

Karen M. Wai, MD

Massachusetts Eye and Ear

11/5/2020



FINANCIAL DISCLOSURES:

None

Ophthalmology Times Research Scholar Honoree Program

MY ROLE IN THIS RESEARCH:

Please answer which of the following portions of the research you participated in:

Conception and design of the work/project
Acquisition of data
Analysis and interpretation of data
Creation and/or critical review of the presentation

Ophthalmology Times Research Scholar Honoree Program

Background: "Why can't I see?"

- Visual acuity is the most widely used test to assess visual function, but may not always accurately describe a patient's self-assessment of visual limitations, especially in those with macular disease
- Contrast sensitivity better correlates with day-to-day activities
- Contrast sensitivity function (CSF) is perhaps a more accurate measurement of visual function but how can we quantify it?



Background: The search for a perfect test

- Absence of a reliable and clinically applicable test has hindered adoption of CSF testing in routine clinical practice
 - Pelli-Robson chart: tests CS only in a single spatial frequency
 - Conventional lab CSF testing: too time consuming
 - Vistech test or functional acuity contrast test: poor test-retest reliability



Ophthalmology Times Research Scholar Honoree Program


Background: How can we apply an intelligent tool for vision testing?

Ophthalmology Times

Research Scholar

The quick contrast sensitivity (qCSF) method

- Use of an adaptive algorithm to select optimal optotypes from >2400 candidate test items
- High sensitivity
- Time efficient
- Maximal information gain



Objective

We evaluate qCSF method in assessing visual function deficits in eyes with macular disease that retained good measured letter acuity compared to healthy control eyes



Methods

- Prospective cross-sectional study
- Patients recruited from December 2016 to February 2020
- Inclusion criteria:
 - Eyes with maculopathy [retinal vein occlusion (RVO), macula-off retinal detachment (RD), dry age-related macular degeneration (AMD), wet AMD] OR healthy control eyes without evidence of maculopathy
 - Retention of good measured visual acuity (VA≧20/30; 0.176 logMAR)
- Exclusion criteria:
 - Visually significant cataracts (>2+ nuclear sclerosis)
 - Presence of any ocular disease (ex. glaucoma)



Methods cont.


- qCSF testing
 - Manifold Contrast Vision Meter (Adaptive Sensory Technology (AST), San Diego, CA)
 - Intelligently selects from 128 possible contrasts, 19 possible spatial frequencies



Sentio

Sentio Response Tablet





Main outcomes:

- Area under the log contrast sensitivity function curve (AULCSF)
- Contrast sensitivity thresholds at six spatial frequencies from 1 – 18 cycles per degrees (CPD)


Results

- 151 eyes from 120 patients with a maculopathy compared to 93 control eyes from 82 patients
 - 27 eyes with RVO from 25 patients
 - 22 eyes with history of macula-off RD from 22 patients
 - 84 eyes with dry AMD from 58 patients
 - 18 eyes with wet AMD from 15 patients




Results

- Presence of maculopathy was associated with:
 - Reduction in AULCSF relative to control eyes
 - Reduction in CS thresholds among all spatial frequencies except at 18 CPD relative to control eyes



Results

- Notably, reductions in CS thresholds were most evident at intermediate spatial frequencies
- Subgroup analysis for very good visual acuity patients (VA≧20/20⁻
 ¹) – similar results



Conclusions: A better assessment of how well our patients can see

- Eyes with macular disease, <u>despite high letter acuity</u>, had significant reductions in contrast sensitivity function relative to healthy controlled eyes
- qCSF is a promising visual function endpoint given its ability to recognize and assess functional deficits in eyes with macular disease that otherwise can go unrecognized by current testing methods



Thank you!

- Dr. John B. Miller
- Harvard Retinal Imaging Lab
- Ophthalmology Times Research Scholar Honoree Program

